


PRODUCT OVERVIEW

Convenient, Fast and Easy-to-Use

# Safe-T-Vue®



## Nonreversible Temperature Indicators for Blood Components



- For safe transport and temporary storage of blood
- Assures temperature compliance and improves quality assurance
- Specifically indicates blood may be reissued when returned\*
- Promotes requirements of accountability during transport
- Proven to improve compliance and preserve valuable blood supplies



**Safe-T-Vue® 6**  
Temperature Indication: 6°C

Use for:

- Refrigerator storage
- Temporary container storage
- For products that may be returned for reissue

Components:

- RBCs
- Whole Blood
- Plasma



**Safe-T-Vue® 8**  
Temperature Indication: 8°C

- Transport in temporary storage containers
- Pneumatic tube transport

- Vaccines
- Biologicals
- Specimens



**Safe-T-Vue® 10**  
Temperature Indication: 10°C

- Transport containers
- Pneumatic tube Transport
- For products that will not be returned or reissued

- RBCs
- Whole Blood
- Plasma
- Biologicals

\*AABB Standards for Blood Banks and Transfusion Services and 21 CFR600.15 Code of Federal Regulations for Temperatures During Shipment for Biological Products



PRODUCT OVERVIEW

Convenient, Fast and Easy-to-Use

# TRANS-VUE™

Nonreversible Temperature Indicators for Transport and Storage Containers

- Monitors temperature for transport and storage in a controlled temperature environment
- Protects shipment integrity of sensitive items for blood banks and laboratories
- Unique, reusable design maintains constant contact with transport environment
- May be reused for up to three months



**TRANS-VUE™ 6**  
Temperature Indication: 6°C

Use for:

- Refrigerator storage
- Temporary container storage
- For products that may be returned for reissue

Components:


- Blood Products
- Biologicals
- Specimens



**TRANS-VUE™ 8**  
Temperature Indication: 8°C

- Transport in temporary storage containers
- Pneumatic tube transport

- Vaccines
- Biologicals
- Specimens



**TRANS-VUE™ 10**  
Temperature Indication: 10°C

- Transport containers
- Pneumatic tube transport
- For products that will not be returned or reissued

- RBCs
- Whole Blood
- Plasma
- Biologicals

An Effective Way to Control Painting of Newborn Umbilical Stumps

## Triple-Dye™

Disposable Pin-Point Applicator with Fiber-Free Bio-Filter Tip™



William Laboratories, Inc.

5 Anngina Drive, Unit B  
Enfield, CT 06082 USA

Toll Free: 800-767-7643

860-749-1350

Fax: 860-749-1351

info@williamlabs.com

